

MAGNIFICO™ Loupes

UNICO® is proud to introduce the MAGNIFICO™ line of magnifying Loupes.

Crisp, clear imaging is provided by these very useful loupes. By enhancing your ability to see the smallest details; these loupes help make sure you can accurately make that critical diagnosis, or get the procedure done correctly.

These attractive glasses feature swing-out loupes and eye protection. When you don't need them for a moment, simply swing them up, and out of your line of sight, then back down again the instant you need them.

The loupes adjust interpupillary distance with one simple twist of a knob.



Working Distance: E2500 & E3500 = 420mm (16.5") +/- 30mm

E2500 Field of View: 100mm – 140mm (2.5" – 5.5")

E3500 Field of View: 60mm – 115mm (2.3" – 4.5")



Each package contains a set of loupes mounted on safety glasses, a lens cloth for keeping them clean, a screw driver for making tension adjustments, lanyard so you can hang them around your neck, and a sturdy, protective case for storage when not in use.

Item #	Description	UOM
E2500	2.5x Loupes, complete with case, lens cloth, screw driver and lanyard, 420 mm (16.5") working distance	Each
E3500	3.5x Loupes, complete with case, lens cloth, screw driver and lanyard, 420 mm (16.5") working distance	Each

UNICO®

UNITED PRODUCTS & INSTRUMENTS

182 Ridge Road, Suite E Dayton, New Jersey 08810

Tel: 800-588-9776 * Fax: 732-274-1151

E-mail: info@unicosci.com * Web: www.unicosci.com

© 2013 UNICO®



MAGNIFICO™ LED Package

Portable LED light option from UNICO offers a lot of light in a very small package.

Offering up to 15 hours of use with a single charge this innovative light source will take you through your work day without slowing you down.

Completely portable with the external long life LIL battery offers 3 different intensity levels.

Pure white light allows user to see everything in true colors.



Color Temperature - 6000-7000K	Brightness - 40000Lux
Beam Color—Pure White	Bulb Life — 50,000 Hours



Each unit is complete with :

- 1 LED headlight with loupes mount
- 1 Battery pack with integrated belt clip
- Battery charger
- Storage case, lens cleaning cloth and user manual.

Item #	Description	UOM
E5000	LED Clip-on headlight, with battery, charger and storage case	each

UNICO®

UNITED PRODUCTS & INSTRUMENTS

182 Ridge Road, Suite E Dayton, New Jersey 08810

Tel: 800-588-9776 * Fax: 732-274-1151

E-mail: info@unicosci.com * Web: www.unicosci.com

© 2013 UNICO®

