

# Porta-Spin™ Mobile Centrifuge from UNICO

- Designed for the professionals on the go
- Plugs into a car 12V power outlet
- 6 or 8 place rotor
- 500-3800RPM variable speed
- 30 minute timer
- Dual lid latch for added safety
- Convenient stainless steel carrying handle
- Optional AC power adapter available
- Optional Tube Holdster™ - side mounted balance tube rack



C829H with Tube Holdster™



## UNICO Professional Digital Non-Contact Photo-Tachometer

- Highly accurate, scientific tool
- Easy to operate
- Standard adhesive reflective tape (included)
- Aim and measure--RPMs instantly display on large LCD screen.
- Includes custom protective cushioned storage case
- Order pn# C800-75

UNICO Model	Item #	Tube Holdster™	Speed	Max. RCF (g)	Angled Rotor
Porta-Spin™	C827	No	500-3800 RPM	1647	6 place
	C829	No			8 place
	C827H	Yes			6 place
	C829H	Yes			8 place



United Products & Instruments, Inc  
[www.unicosci.com](http://www.unicosci.com) \* [info@unicosci.com](mailto:info@unicosci.com)

182 Ridge Road, Suite E Dayton, NJ 08810  
 Tel: 732-274-1155 \* Fax: 732-274-1151

Porta-Spin™, Tube Holdster™ are Registered Trademarks of United Products & Instruments, Inc.