



- Provides stabilization, cold therapy, and compression
- Helps relieve pain and swelling after an injury
- Contoured design fits comfortably on left or right ankle
- Combination gel/air bladders incorporate cold therapy and proper support
- Double sealed bladders prevent leakage

Sized by US shoe size

MODEL #	SIZE	MEN'S	WOMEN'S
AKL-6371-WH-SML	S	3-5	4-6
AKL-6371-WH-MED	M	6-7	7-8
AKL-6371-WH-LRG	L	8+	9+



Help Relieve Pain & Swelling

The Swede-O™ Versi-Splint™ provides a great combination of medial and lateral ankle stabilization, cold therapy, and compression to help relieve pain and swelling following acute ankle injuries.

Upper bladder contains air while lower bladder contains gels. This combination provides proper compression with cold therapy.

Lightweight & Durable

These durable, vinyl bladders are double sealed to prevent leakage. The ultralight, contoured rigid shell fits comfortably over the right or left ankle. Thin heel pad is comfortable and easily fits in most shoes. To utilize the splint for cold therapy, simply place the product in the freezer for 30 minutes prior to use.

Specifications

Product Selection:

AKL-6371 Swede-O™ Versi-Splint™

Product Weight: .6 lb

Dimensions: 9" x 11.5" x 3"

Case Weight: 8 lbs

Case Dimensions: 11" x 14" x 13"

Case Quantity: 12

Available From:

