

800-525-9339 • 651-674-8301  
www.swedeo.com

- Easy to apply and adjust
- Removable side stabilizers for additional support and protection
- Full elastic back reduces irritation and potential for blistering
- Bilateral design fits both right and left ankle
- Industry leading 12 month warranty, satisfaction guarantee

Sized by US shoe size

BLACK MODEL #	WHITE MODEL #	SIZE	MEN'S	WOMEN'S
AKL-6332-BK-SML	AKL-6332-WH-SML	SM	6 - 7 1/2	7-9 1/2
AKL-6332-BK-MED	AKL-6332-WH-MED	MED	8-10 1/2	10-11 1/2
AKL-6332-BK-LRG	AKL-6332-WH-LRG	LRG	11-12 1/2	12-13 1/2
AKL-6332-BK-1XL	AKL-6332-WH-1XL	XL	13-14 1/2	14-15 1/2



## Comfortable & Easy to Use

The Swede-O Easy Lok™ is simple to apply and adjust using the two strap closure system: Ideal for people with arthritic hands and those who have difficulty applying a lace-up style brace.

The Easy Lok's removable side stabilizers offer additional support for weak or injured ankles. The full elastic back is comfortable against the Achilles tendon and helps reduce irritation and potential for blistering. The durable vinyl construction and industry leading 1 year warranty provide peace of mind that you're giving your ankle the best protection. This brace fits comfortably in most shoes.

## Indications for Use

Use the Easy Lok after injury or to protect your ankle from future injury.

This ankle brace helps prevent and relieve aches, pains, sprains, and other ankle conditions.



Available From:

## Specifications

### Product Selection:

#AKL-6332 Swede-O Easy Lok™

**Product Weight:** .35 lb

**Dimensions:** 9" x 11.5" x 3"

**Case Weight:** 7 lbs

**Case Dimensions:** 11" x 14" x 13"

**Case Quantity:** 12

**Materials:** Vinyl, polyester, nylon, polyurethane foam, HDPE

