

- Easily control support with external figure eight straps
- Durable, lightweight ballistic nylon is ideal for sports
- Cool, comfortable, breathable tongue
- Full elastic back reduces irritation and potential for blistering
- Bilateral design fits right or left ankle
- Industry leading 1 year warranty, satisfaction guarantee
- Helps lock the ankle in a protected position

Sized by US shoe size

MODEL #	SIZE	MEN'S	WOMEN'S
AKL-6338-BK-1XS	XS	3-5 1/2	4-6 1/2
AKL-6338-BK-SML	SM	6 - 7 1/2	7-9 1/2
AKL-6338-BK-MED	MED	8-10 1/2	10-11 1/2
AKL-6338-BK-LRG	LRG	11-12 1/2	12-13 1/2
ALK-6338-BK-1XL	XL	13-14 1/2	14-15 1/2

Simulates Ankle Taping

The Swede-O Multi-Sport™ ankle brace has external figure eight straps which allow you to easily control and customize the amount of support, as well as simulate professional taping.

The top strap covers the laces and locks the figure eight straps in position. The full elastic back is comfortable against the Achilles tendon and helps reduce irritation and the potential for blistering. The breathable tongue helps keep the foot cool and comfortable.

The Multi-Sport is backed by an industry leading 1 year warranty and constructed from a durable ballistic nylon which is comfortable thin, and lightweight.



Available From:



LATEX FREE

SWEDE-O
by Core Products
INTERNATIONAL, INC.

Specifications

Product Selection:

#AKL-6338 Swede-O Multi-Sport™

Product Weight: .35 lb

Dimensions: 9" x 11.5" x 3"

Case Weight: 7 lbs

Case Dimensions: 11" x 14" x 13"

Case Quantity: 12

Materials: Polyester, nylon, polyurethane foam, neoprene, brass

