

- Quick and easy speed lacing
- Integrated stabilizer conforms to the ankle for superior support and protection
- Low profile - Fits in most any shoe
- Support of a rigid brace, comfort of a lace up brace
- Industry leading 12 month warranty, satisfaction guarantee

Sized by US shoe size

MODEL #	SIZE	MEN'S	WOMEN'S
AKL-6337-BK-1XS	XS	3-5 1/2	4-6 1/2
AKL-6337-BK-SML	SM	6-7 1/2	7-9 1/2
AKL-6337-BK-MED	MED	8-10 1/2	10-11 1/2
AKL-6337-BK-LRG	LRG	11-12 1/2	12-13 1/2
AKL-6337-BK-1XL	XL	13-14 1/2	14-15 1/2



Support of a Rigid... Comfort of a Lace-up

The Swede-O Tarsal Lok® features a speed lacing system for quick and easy application. The integrated design uses body heat to conform the stabilizer panel to the shape of your ankle and foot, providing a custom fit for superior protection and support.

The breathable tongue helps keep the foot cool and comfortable. The Tarsal Lok is backed by an industry leading 1 year warranty and constructed from a durable, abrasion resistant material.

Indications for Use

Use the Tarsal Lok after injury or to help protect your ankle from future injury. This ankle brace helps prevent and relieve aches, pains, sprains, and other ankle conditions.



Available From:



Specifications

Product Selection:

#AKL-6337 Swede-O Tarsal Lok®

Product Weight: .35 lb

Dimensions: 9" x 11" x 3"

Case Weight: 7 lbs

Case Dimensions: 11" x 14" x 13"

Case Quantity: 12

Materials: Polyester, nylon, polyurethane foam, neoprene, brass, polypropylene

Patent: 5,741,222

