

800-365-3047 • 715-294-2050
www.coreproducts.com

- Deep, penetrating massage with a twist of the switch
- Helps reduce pain, muscle tension, and stiffness
- Variable speeds provide soothing massage or more invigorating massage
- An industry standard among professional massage therapists
- Helps stimulate the circulatory system
- Fosters healing of strained muscles and sprained ligaments



Therapeutic Relief

Making bodies feel better since 1966, the original Jeanie Rub® delivers invigorating 1100-1400 rpm massaging action with a twist of the switch.

The orbital action provides a softer and deeper massage than “percussion” type massagers. Convenient fingertip control allows you to easily adjust speed – lower speeds for soothing relief and faster speeds for a more invigorating experience. Enhance your massage techniques while reducing stress on your wrists and hands.

The Jeanie Rub is ideal for personal health, fitness and sports massage. It helps to naturally manage pain associated with chronic or overuse injuries by easing tension, promoting flexibility and helping to improve circulation. As an addition to any

sports training regimen, it helps to decrease the risk of injury and reduce workout recovery time.

High Quality Construction

The Jeanie Rub is designed to withstand years of use. All Jeanie Rub massagers are built in the USA to ANSI Standards and are CSA certified, ensuring they have been evaluated to meet or exceed consumer safety standards. The cast aluminum construction with red and black powder coat finish looks great and resists abrasions. Pad measures 4” x 9” (10cm x 23cm). Jeanie Rub accessories sold separately.

Available From:



Specifications

Product Selection:

PRO-3401 Jeanie Rub Massager®

Product Weight: 7 lb

Dimensions: 9.5” x 7.75” x 4.5”

Case Weight: 28 lbs

Case Dimensions: 19” x 11” x 8.5”

Case Quantity: 4

800-365-3047 • 715-294-2050
www.coreproducts.com



Paraspinal Accessory
ACC-896



Extremity Accessory
ACC-897



Wall Hanger
ACC-884

Professional Package Includes:

- Jeanie Rub Variable Speed Massager
- Paraspinal Accessory
- Extremity Accessory
- Carrying Bag
- Fleece Cover

Professional Package
PRO-3405



Terry Cover
ACC-893



Sure-Chek Cover
ACC-892



Fleece Cover
ACC-882



Sheepskin Cover
ACC-886



Carrying Bag
ACC-885

Specifications

Product Selections:

#PRO-3405: Professional Package
#ACC-896: Paraspinal Accessory
#ACC-897: Extremity Accessory
#ACC-884: Wall Hanger
#ACC-893: Terry Cover (3 pk)
#ACC-892: Sure-Chek Cover (3 pk)
#ACC-882: Fleece Cover
#ACC-886: Sheepskin Cover
#ACC-885: Carrying Bag

Dimensions:

10.25" x 13.5" x 12.5"
10.25" x 5" x 2"
10.25" x 5" x 4.5"
14" x 5" x 5"
10" x 7"
9" x 4"
9.5" x 5.5"
10" x 8"
12.5" x 8.5" x 8.5"

Wt:

11 lb
.15 lb
.2 lb
1.5 lb
.2 lb
.15 lb
.1 lb
.3 lb
1 lb

Case Qty:

1
6
6
1
3
3
12
6
6

Case Wt:

11 lbs
2 lbs
2 lbs
3 lbs
2 lbs
2 lbs
4 lbs
3 lbs
8 lbs



Available From:

