

T340 Colored Closures

Six styles of caps to choose from and two sealing types: washer seal and lip seal.

The cap is molded with a deep internal lip that fits snugly against the interior wall of the tube thus preventing the contents from coming in contact with the seal or threads, thereby reducing the chances of sample contamination. The cap's high profile facilitates manipulation especially in aseptic procedures and can remain attached (T332 & T336 Series) to the tube in order to eliminate mix-ups and contamination. All tubes and closures are manufactured in a clean environment.

Closures can be COLOR CODED by the use of T345 Series Colored Capinsert™ inserted on the top of the closure. This is accomplished without removing the cap. Colored caps are also available in all models as listed below.



With WASHER SEAL and without loop. Made of polypropylene.

LIP SEAL and without loop. Made of polyethylene.



Cat. #	Cat. #	Color	Qty/Pk
T340 NOS	T340 NLS	Natural	1000
T340 BOS	T340 BLS	Blue	1000
T340 GOS	T340 GLS	Green	1000
T340 LOS	T340 LLS	Lilac	1000
T340 OOS	T340 OLS	Orange	1000
T340 ROS	T340 RLS	Red	1000
T340 YOS	T340 YLS	Yellow	1000
T340 WOS	T340 WLS	White	1000
T340 BROS	T340 BRLS	Brown	1000



With WASHER SEAL and loop. Made of polypropylene.

LIP SEAL and loop. Made of polyethylene.



Cat. #	Cat. #	Color	Qty/Pk
T340 NOSL	T340 NLSL	Natural	1000
T340 BOSL	T340 BLSL	Blue	1000
T340 GOSL	T340 GLSL	Green	1000
T340 LOSL	T340 LLSL	Lilac	1000
T340 OOSL	T340 OLSL	Orange	1000
T340 ROSL	T340 RLSL	Red	1000
T340 YOSL	T340 YLSL	Yellow	1000
T340 WOSL	T340 WLSL	White	1000
T340 BROS	T340 BRLS	Brown	1000

The following closures have a flat top to accommodate automatic capping machines in packaging industries.



With WASHER SEAL and without loop. Made of polypropylene.

LIP SEAL and without loop. Made of polyethylene.



Cat. #	Cat. #	Color	Qty/Pk
T340 NOSFT	T340 NLSFT	Natural	1000
T340 BOSFT	T340 BLSFT	Blue	1000
T340 GOSFT	T340 GLSFT	Green	1000
T340 LOSFT	T340 LLSFT	Lilac	1000
T340 OOSFT	T340 OLSFT	Orange	1000
T340 ROSFT	T340 RLSFT	Red	1000
T340 YOSFT	T340 YLSFT	Yellow	1000
T340 WOSFT	T340 WLSFT	White	1000
T340 BROSFT	T340 BRLSFT	Brown	1000

T345

Color Coding CAPINSERT™



Made of polypropylene

The Capinsert™ is used to color code a Microwtube® and a multitude of other Simport products according to your specific needs. It is inserted on top of the closure and has a write-on frosted area for sample identification.

Cat. #	Color	Cat. #	Color	Qty/Bag
T345B	Blue	T345P	Pink	500
T345GY	Gray	T345R	Red	500
T345G	Green	T345V	Violet	500
T345L	Lilac	T345W	White	500
T345O	Orange	T345Y	Yellow	500
* Blue, lilac, red, yellow and white		T345AS	Assorted*	500