

**Tamper Evident**



## M960 HISTOTAINER™ I

### TAMPER EVIDENT Prefilled Specimen Containers, 50% Filled with 10% Neutral Buffered Formalin

Container made of polypropylene  
Closure made of polyethylene

Especially designed for collection, transport and storage of histology specimens, Simport® offers shatter resistant polypropylene containers, eliminating most problems of leakage and evaporation. The magenta lids are ribbed for easy opening when hands are wet or gloved while the jars are stackable for easy, safe storage and translucent to allow specimens to be viewed without opening. The closures are manufactured from virgin polyethylene with a unique integrated leak-resistant seal. Container's vertical walls offer excellent rigidity.

The uniqueness of the HistoTainer™ I is that it incorporates an innovative

tamper evident screw cap ensuring your peace of mind during transport or storage situations where someone might have manipulated the specimen without your prior knowledge. Can also be used without using the tamper evident locking mechanism. For color coding purposes, use a Capinsert™ (see T345 Series) on top of closure. Ten different colors are available.

Both containers and caps are manufactured without the use of plasticisers or mold release agents. All material used in manufacture are free from latex. All containers are 95 kPa compliant. Available in many sizes from 20 to 120 ml. Packaging in trays of 24, cases of 96.

The Simport® HistoTainer™ I is half filled with 10% Neutral Buffered Formalin as a fixative. 10% Neutral Buffered Formalin penetrates quickly, but fixes slowly. The Simport® Formalin is enhanced by a buffering capacity optimizing histological results by light microscopy and immunohistochemistry.

Cat. #	Packaging	Volume	Qty/Tray	Qty/Cs
M960-20FMA	Tamper Evident – Internal Trays	20 ml	24	96
M960-40FMA	Tamper Evident – Internal Trays	40 ml	24	96
M960-60FMA	Tamper Evident – Internal Trays	60 ml	24	96
M960-90FMA	Tamper Evident – Internal Trays	90 ml	24	96
M960-120FMA	Tamper Evident – Internal Trays	120 ml	24	96



## How to use the HistoTainer™ I



Remove Screw Cap.



Place sample in container.



Push up the attached tab on side of container. Make sure it firmly clicks in place.



Screw cap completely on container.



When opening the container, the tamper evident ring will detach itself from the cap.



The HistoTainer™ can also be used without the tamper evident feature.



As a cassette holder, the 40 ml HistoTainer™ M960-40FMA containing 20 ml of 10% Formalin, is the ideal size container to transport up to four tissue samples pre-inserted in processing / embedding cassettes.



M960 Series packaged in sturdy cardboard boxes with handles for easy carrying.

## Anatomy of the HistoTainer™ I

1. Ridges around base offer a better grip during opening and closing.
2. Warning label has space for patient identification.
3. Insertion of a Capinsert™ allows color coding identification of contents.
4. Molded ridges around lid make it easy to open and close.
5. Tamper evident sealing ring for better sample protection.
6. Specially designed locking tab to ensure a perfect tamper evident seal.
7. 10% Neutral Buffered Formalin helps to protect sample integrity.

