



**C577**

## SecurTainer™ III Specimen Containers

Container made of polypropylene / Closure made of polyethylene



For IVD use

Chemically resistant and shatterproof, The SecurTainer™ III Specimen Containers are supplied with a leakproof screw cap, particularly important when transporting hazardous material. Ideal for transport and storage of urine, sputum and most liquids or particulate samples.

Containers are uniquely stackable, shatter resistant and are manufactured from virgin, high-clarity polypropylene. Available in many sizes from 20 to 120 ml. Containers and caps packaged separately in bags of 100. Non sterile. Functional temperature range: -90 °C to +100 °C.

### Anatomy of a SecurTainer™ III

1. Ridges around base offer a better grip during opening and closing.
2. Insertion of a Capinsert™ allows color coded identification of contents.
3. Molded ridges around lid make it easy to open and close.



Cat. #	Volume	Qty/Pk	Qty/Cs
C577-20W	20 ml	100	500
C577-40W	40 ml	100	500
C577-60W	60 ml	100	500
C577-90W	90 ml	100	400
C577-120W	120 ml	100	300

**T345**

## Color Coding CAPINSERT™



Made of polypropylene

The Capinsert™ is used to color code a multitude of Simport® products according to your specific needs. It is inserted on top of the closure and has a write-on frosted area for sample identification.

Cat. #	Color	Cat. #	Color	Qty/Bag
T345B	Blue	T345P	Pink	500
T345GY	Gray	T345R	Red	500
T345G	Green	T345V	Violet	500
T345L	Lilac	T345W	White	500
T345O	Orange	T345Y	Yellow	500
		T345AS	Assorted*	500

\* Blue, lilac, red, yellow and white