

URISAFE®

24-HR URINE COLLECTION CONTAINER

Container made of high-density polyethylene
Cap made of polypropylene

Uniquely designed to be
the most user friendly



For IVD use



On the 4 liter model, urine can be poured into a tube by simply tilting the container forward without having to lift it.

- 1 The unique snap valve pour spout is easily popped open and offers dripless pouring
- 2 Specially designed leakproof screw cap with liner for safer transport
- 3 Large centrally located handle for good balance and ease of manipulation
- 4 Metal free and resistant to hydrochloric acid
- 5 Sample volumes can be read in 100 ml increments

The innovative cap:

- Snap valve pouring spout incorporated
- Leakproof: a cap liner ensures safe sample transport

The unique snap valve pour spout:

- Easily popped open or pushed shut
- Offers dripless pouring
- Controlled flow rate for better handling of poured volumes
- Reduced risks of aerosol contamination when pouring
- Eliminates splashing and exposure to hazardous body fluids

Low form design:

- Convenient refrigerator storage with less wasted space
- Anatomically designed for ease of patient use

Large central handle:

- Can be gripped comfortably with three fingers
- Permits pouring off of samples with ease and reduces fatigue

Chemically resistant:

- Metal, latex, zinc and fluorescence free
- Resistant to hydrochloric acid.
- Can be subjected to freezing, thawing and EtO gas sterilisation without causing changes in materials or physical appearance

Graduated vertically and horizontally:

- Graduations are easy to read
- Place container upright and read in 100 ml increments
- Can be gamma sterilized

The URISAFE® Urine Collection Containers are USER FRIENDLY to both patients and laboratory personnel. Two sizes are available: 3 and 4 liters.

Cat. #	Vol. (Liters)	Dimensions (cm)	Qty/Cs
B350	3	11.5 x 24.5 x 16.0 H	40
B350-4L	4	11.5 x 24.5 x 20.1 H	30



B350-4HCL Urisafe with HCl Preservative

This model offers the convenience of being pre-filled with 30 ml 6N HCl. All patient instructions are clearly detailed on label.

Each container is individually wrapped.

Cat. #	Vol. (Liters)	Dimensions (cm)	Qty/Cs
B350-4HCL	4	11.5 x 24.5 x 20.1 H	30