

# AirLife® soft cushion cannulas

AirLife® soft cushion cannulas include lariat tubing that is approximately 37 percent\* softer than our standard AirLife cannulas. This soft tubing helps alleviate chafing behind the ears, reducing skin irritation.

## Features:

- Optimum drapability to secure placement
- Soft lariat to reduce need for foam ear covers
- Wide variety of lengths
- Soft design to optimize patient comfort
- Crush-resistant tubing to prevent occlusion and kinking
- Latex-free materials
- Not made with DEHP

Catalog no.	Description	Qty. (cs.)
SFT2611	Neonate, with 7' O <sub>2</sub> tubing	50
SFT2601	Infant, with 7' O <sub>2</sub> tubing	50
SFT2602	Pediatric, with 7' O <sub>2</sub> tubing	50
SFT2692	Pediatric, with 7' O <sub>2</sub> U/Connect-It® tubing	50
SFT2604	Adult, with 4' O <sub>2</sub> tubing	50
SFT2600	Adult, with 7' O <sub>2</sub> tubing	50
SFT2614	Adult, with 14' O <sub>2</sub> tubing	50
SFT2625	Adult, with 25' O <sub>2</sub> tubing	25
SFT2699	Adult, with 7' O <sub>2</sub> U/Connect-It tubing	50
SFT2614U	Adult, with 14' O <sub>2</sub> U/Connect-It tubing	50

 **WARNING**—U.S. Federal Law restricts this device to sale by or on the order of a physician.

CareFusion  
Vernon Hills, IL

[carefusion.com](http://carefusion.com)

\*Based on a material durometer difference of 30 using the Shore A hardness scale.

© 2015 CareFusion Corporation or one of its affiliates. All rights reserved. AirLife, UConnect-It, CareFusion and the CareFusion logo are trademarks or registered trademarks of CareFusion Corporation or one of its affiliates AL4926 (0415/300)

