

MicroGard™ II Pulmonary Function Filter

The MicroGard II pulmonary function filter acts as an effective barrier for the protection of patients, staff and instrumentation from microbial contamination.

MicroGard II stands for patient and operator safety without compromising on system performance and clinical output.

Key features:

- 99,999% viral and bacterial efficiency against cross contamination
- Low resistance to airflow
- Minimal dead space 55 mL
- Filter housing and filter fleece tested for biocompatibility
- 2 filters and 2 filter kits configuration versions
- Validated for compliance with the complete CareFusion Respiratory Diagnostics product line



CareFusion

The MicroGard II filter series

- **MicroGard IIB** with integrated oval shaped mouthpiece suitable for adult and pediatric use.
- **MicroGard IIC** with round shaped connector for use with special lung function mouthpieces, ideal for gas dilution and bodyplethysmographic measurements.

After extensive testing at independent labs such as Nelson Laboratory the MicroGard II bacterial/viral filter confirms to offer the right balance between optimized filtration efficiency and minimal resistance (Delta P) to airflow. Based on strict CareFusion verification and validation processes, the MicroGard II does not affect equipment measurement characteristics while protecting patient, operator and equipment from cross-contamination. Therefore MicroGard II filters are the filters of choice for use with CareFusion JAEGER®, Vmax® and Micro Medical PFT instrumentation.

Specifications:

- **Inspiratory resistance:** < 0.04 kPa/(L/s) at 1 L/s (< 0.4 cmH₂O/(L/s) at 1 L/s)
- **Expiratory resistance:** < 0.04 kPa/(L/s) at 1 L/s (< 0.4 cmH₂O/(L/s) at 1 L/s)
- **Filter Efficiency against cross contamination:**
 - Bacterial > 99.999% (based on Nelson Test Report 10003754)
 - Viral > 99.999% (based on Nelson Test Report 10003754)

Filter kit MG IIB
(MicroGard IIB filter + noseclip)




Filter kit MG IIC
(MicroGard IIC filter
+ disposable rubber mouthpiece + noseclip)



Ideal for gas dilution testing
and for poorly cooperating subjects

Ordering information, P/N and description:

- V-892381:** MicroGard IIB (single filter, integrated oval shaped mouthpiece) - Box of 50
- V-892384:** MicroGard IIC (single filter, round connector) - Box of 50
- V-892391:** Filter kit (MicroGard IIB + noseclip) - Box of 80
- V-892392:** Filter kit (MicroGard IIC + disposable rubber mouthpiece + noseclip) - Box of 80
- V-861449:** Silicon adapter, oval, for volume calibration using MicroGard IIB
- 852740:** Silicon adapter, round, for volume calibration using MicroGard IIC
- 36-MFA1010:** adapter MicroGard to Micro Medical spirometers, pack of 10
- 36-MFA1050:** adapter MicroGard to Micro Medical spirometers, pack of 50

 CareFusion Germany 234 GmbH
Leibnizstrasse 7
97204 Hoechberg
Germany

+49 931 4972-0 tel
+49 931 4972-423 fax



CareFusion
22745 Savi Ranch Parkway
Yorba Linda, CA 92887

800.231.2466 toll-free
714.283.2228 tel
714.283.8493 fax

U.K. Sales
CareFusion UK 236 Ltd
The Crescent, Jays Close
Basingstoke, RG22 4BS, UK

+44 (0) 1256 388599 tel
+44 (0) 1256 330860 fax

carefusion.co.uk/microgard

© 2016 CareFusion Corporation or one of its affiliates.
All rights reserved. CareFusion, CareFusion Logo, MicroGard, JAEGER and Vmax are
trademarks or registered trademarks of CareFusion Corporation or one of its affiliates.
All trademarks are property of their respective owners. V-791718 Rev.4

