

# Electrotherapy

## InTENSity™ 10 Digital TENS Unit

The InTENSity™ 10 Digital TENS device allows patients to select from ten pre-set body icons for targeted therapy, making it the easiest and most effective TENS device available on the market today. With the press of two buttons the easy-to-use device is ready to operate effectively based on the area of the body that is being treated. Comes complete with AC adaptor for effective therapy on the go.



### Product Features

- Simple to operate, portable device for pain management
- 2-button operation
- Device has ten pre-set body part selections making it easy for patients to target certain areas in pain
- Body part selections include: neck, shoulder, mid back, low back, elbow, wrist, hip, knee, ankle and foot
- Individual body part selections have a customized treatment program based on researched protocols
- Comes complete with carrying case, two sets of lead wires, 9-volt battery, AC adaptor, manual and pack of four self-adhesive reusable electrodes
- Power Source: Standard 9-volt battery or AC adaptor
- Dimensions: 4.5" x 2.55" x 0.9"
- Weight: 0.28 lb. with battery
- One year warranty
- HCPCS Code: E0730\* for TENS Unit
- HCPCS Code: A4595\* for electrotherapy supplies



### Product Specifications:

M= Modulate F= Fixed

	Neck	Low Back	Mid Back	Shoulder	Elbow	Wrist	Knee	Hip	Foot	Ankle
Pulse Width	M 100-150	F 260	F 100	F 260	F 100	F 100	M 100-150	F 200	M 100-160	F 100
Pulse Rate	M 60-100	M 50-80	M 100-150	M 80-100	F 100	F 260	F 120	M 100-150	M 60-100	F 100

Item #	Description	Qty
D11010	InTENSity™ 10	each

\* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



**To Order:**  
 (800) 3-ROSCOE  
 fax: (440) 572-4261  
 roscoemedical.com