

Lumin™



The Lumin helps CPAP patients stay healthy and sleep better.

To learn more, contact:



Introducing

Lumin™

CPAP Sanitizing System

The healthy choice in CPAP accessory disinfection with simple one-touch operation.

Lumin is the easiest and fastest way to disinfect a CPAP mask, reservoir and other accessories. The Lumin uses UV-C light to safely sanitize and works with a disinfection cycle time of 5 minutes, no harmful ozone, and 99.9% kill rate for harmful bacteria, viruses, mold and fungus. It is the ideal system for daily sanitizing of your CPAP equipment.

The Lumin uses UV light to safely and naturally sanitize your CPAP mask and accessories.



The Lumin was developed out of the need for a product that would not damage the CPAP device itself during disinfection and had zero risk of exposing a respiratory patient to a lung irritant (no harmful Ozone).

Lumin benefits:



Lumin features:

- Sleek design
- Small footprint
- Safe UV-C light
- Simple, one-button operation
- Status indicator lights
- Automatic shut-off
- 24-month warranty

How the Lumin works:

UV-C light is considered to be the safest disinfection option on the market, as Ozone used in other systems has been shown by the EPA to cause irritation of the lungs and respiratory system at relatively low amounts. The Lumin's UV-C light source with a 99.9% kill rate of bacteria and germs, significantly reducing the risk of infection from an improperly cleaned CPAP.

Simply place your mask, reservoir and accessories in the Lumin drawer and press the on/off button to begin the 5-minute disinfection cycle. Full and complete sanitizing is achieved without the use of messy chemicals or water.