



## TENS 7000™ 2<sup>nd</sup> Edition (OTC)

### TENS THERAPY

The TENS 7000 2<sup>nd</sup> Edition is one of the most popular, and best valued, digital TENS units on the market and is now available without a prescription. This digital TENS unit is proven to be among the most durable and easy-to-use TENS units providing drug-free pain relief for sore or aching muscles of the lower back, shoulder, arm, leg due to strain from exercise and normal activities. The upgraded digital processor is more precise and accurate than ever before. This unit is great for users seeking a digital device that is also equipped with physical amplitude knobs. FDA cleared for over-the-counter sale.



#### PRODUCT FEATURES:

- Providing drug-free pain relief for sore or aching muscles of the lower back, shoulder, arm, leg due to strain from exercise and normal activities
- Modes: Burst, Normal, Modulation, SD1 and SD2
- Dual channel for use with two or four electrodes
- Pulse width and a pulse frequency are adjustable
- System includes device, lead wires, four self-adhesive reusable electrodes, 9-volt battery, hard plastic carrying case, quick start guide and instruction manual
- Compliance meter, 60 records up to 999 hours
- Low battery indicator
- One year warranty
- FDA cleared for over-the-counter sale

#### SPECIFICATIONS:

Channel:	Dual, isolated b/w channels
Power Source:	9-volt battery
Waveform:	Asymmetric Biphasic
Pulse Amplitude:	Adjustable, 0-100 mA at 500 ohm load each channel
Pulse Rate:	Adjustable, 2 to 150 Hz
Pulse Width:	Adjustable, 50~300 $\mu$ s, 10 $\mu$ s/step
Modes:	Burst Normal Modulation Strength Duration - SD1 & SD2
Timer:	5 to 60 minutes with auto shutoff or continuous

Item #	Description	Qty
DT7202	TENS 7000 2nd Edition (OTC)	each
EF2020TC2	Premium Tan Cloth Electrode, Foil Bag	4/pk