

Electrotherapy

TopTENS™ Pain Relief System

The TopTENS™ Pain Relief System is a doctor recommended, over-the-counter (OTC) device that provides prescription-strength pain relief associated with sore and aching muscles in the upper and lower extremities due to strain from exercise or normal household and work activities. This easy-to-use device comes equipped with eight preset programs allowing you to customize your therapy and target your pain. Comes complete with three AAA batteries, lead wires and four reusable electrodes for effective therapy on the go. A great retail item as it is FDA cleared for over-the-counter sale.



Product Features:

- FDA cleared for over-the-counter sale (no prescription required).
- TENS Technology stimulates body's own natural pain killers.
- High intensity output is greater than other portable TENS devices
- Simple-to-operate, portable device for pain management.
- Device has eight pre-set programs allowing you to customize your therapy and target your pain.
- Comes complete with TopTENS™ device, two lead wires, 3 AAA batteries, user manual and pack of four self-adhesive reusable electrodes.



Technical Specifications:

- Channel: Dual, isolated between channels
- Pulse Intensity: Adjustable 0 - 100A peak into 500 Ω load each channel
- Pulse Rate: 2 - 60 Hz
- Pulse Width: 156 - 260 μ s
- Timer: 5-60 min (optional 5 minute increments) and continuous mode
- Wave Form: Symmetrical Bi-phasic square pulse
- Power Source: 3 AAA batteries / 4.5 volt batteries
- Dimensions: 4.5" x 2.25" x 0.9"
- Weight: 0.25 lb. with battery
- Warranty: One year warranty from date of purchase

Item #	Description	Qty
DT6030	TopTENS™ Pain Relief System	each
EIP2020BC2	2" x 2" Reusable Electrodes, Square Blue Cloth, Poly Bag	4/pk
WW6000	Lead Wires For Use With TopTENS™ Pain Relief System	2/pk

*This device is FDA cleared for over-the-counter sale for pain relief associated with sore and aching muscles in the upper and lower (ankles, calves, knees, thigh, hip, hand, wrist, arm, elbow and shoulder) extremities.



To Order:

(800) 3-ROSCOE
fax: (440) 572-4261
roscoemedical.com