

Professional Devices

US 1000™ 3rd Edition Portable Ultrasound Unit

The US 1000™ offers effective ultrasound therapy by generating deep ultrasonic waves within body tissues to treat pain, muscle spasms and joint stiffness. The portable unit comes as a complete kit with carrying case, AC adaptor, ultrasound gel and manual.



Product Features:

- Therapeutic ultrasound device that generates deep heat to effected area
- Device must be used with ultrasound gel to conduct the ultrasound energy to the skin
- Three power output settings
- Comes complete with ultrasound, AC adaptor, carrying case, ultrasound gel and instruction manual
- Six month warranty on unit

Product Specifications:

- Power Source: AC adaptor
- Output Frequency: 1.0 MHz +/- 10%
- Temporal Maximum Power: 4 W +/- 20%
- Output Intensity: Low 1.6 W, Medium 2.6 W, High 4.6 W
- Duty Cycle: Low 8%, Medium 14%, High 29%
- Output Mode: Pulse
- Pulse Width: 8% 05 ms +/- 10%
- Repetition Rate: 150 Hz +/- 10%
- Timer Settings: 30 minutes with auto off

Item #	Description	Qty
DUI025	US 1000™ 3rd Edition Portable Ultrasound Unit	each
LC2828	TheraSonic™ Ultrasound Gel, 2 ounces	each
LS5288	TheraSonic™ Ultrasound Gel, 8 ounces	each
LS5266	TheraSonic™ Ultrasound Gel, 16 ounces	each

* Prescription Required: This is a Class II Medical Device. Federal law restricts these devices to sale by or on the order of a licensed health care professional.

* Providers are responsible for determining appropriate billing codes for Medicare Program claims. This information is included for your convenience and is not intended, nor should it be considered billing advice.



To Order:

(800) 3-ROSCOE
fax: (440) 572-4261
roscoemedical.com