



MATERIAL SAFETY DATA SHEET

Product name: Ultrasound Gel

Model Number: (Ultrasound gel comes with in following devices)

- DU1025
- DU3035
- DU6012
- DU3038
- DU3037

Date prepared: October 28, 2014

1. COMPANY IDENTIFICATION:

Manufacturer: Roscoe Medical Inc.
21973 Commerce Pkwy
Strongsville, OH. 44149

Department: Quality Department

Telephone Number: 800-633-4839

Fax Number: 650-697-3596

2. HEALTH HAZARD INFORMATION

This gel is not dangerous.

Potential Health Effects:

Eye: Transient irritation

Skin: Redness, itching, dryness

Inhalation: None expected

Ingestion: None expected

3. COMPOSITION / INFORMATION ON INGREDIENTS

Main ingredients: Proprietary

4. FIRST-AID MEASURES

Eye Contact:

Flush with water for 10-15 minutes. If necessary, consult medical help.

21973 Commerce Parkway
Strongsville, OH 44149

Toll free: 800-3-ROSCOE (376-7263)
main: 440-572-1962
fax: 440-572-42691

www.roscoemedical.com



MATERIAL SAFETY DATA SHEET

Skin Contact:

Wash with water.

Inhalation:

Not applicable

Ingestion:

Flush mouth with plenty of water. Consult medical help if necessary.

5. FIRE AND EXPLOSION DATA

Flash Point: Not determined

Upper Exposure Limit: None

Lower Exposure Limit: None

Extinguishing Media: All standard firefighting media

Unusual Hazards: N/A

6. ACCIDENTAL RELEASE MEASURES

Clear gel with water and absorbent materials. Cleaning, the product may cause slip hazards.

7. HANDLING AND STORAGE

Storage & Handling Precautions: Store in a cool, dry area. Avoid eye contact or prolonged skin contact.

8. EXPOSURE CONTROLS / PERSONAL PROTECTION

Respiratory Protection: Not required

Eye Protection: Not applicable

Hand Protection: Not required

Skin Protection: Not applicable

Ventilation: Not required

9. PHYSICAL AND CHEMICAL PROPERTIES

PH: 5.5-8.0

Density: 1.07 g/cc

Boiling Range: decomposes at 635°F

Volatile Organic Compounds: None

Solubility in water: 40%

Percent Solids: <5%

Vapor Density: N/A

Vapor Pressure: N/A

Appearance: Clear gel or light color

Odor: Bland odor

21973 Commerce Parkway
Strongsville, OH 44149

Toll free: 800-3-ROSCOE (376-7263)
main: 440-572-1962
fax: 440-572-42691

www.roscoemedical.com



MATERIAL SAFETY DATA SHEET

10. STABILITY AND REACTIVITY

Stability: Stable

Hazardous Decomposition Products: Burning can produce carbon monoxide and/or carbon dioxide and nitrogen oxides.

Hazardous Polymerization: Will not occur

Incompatibility: Strong oxidizing agents

11. TOXICOLOGICAL INFORMATION

Skin irritation: Not determined

Oral toxicity: Not determined

Eye Irritation: Not determined

Known/Suspected Carcinogens: None

12. ECOLOGICAL INFORMATION

N/A

13. DISPOSAL CONSIDERATIONS

Follow applicable federal, state and local regulations. Do not allow to contaminate ponds, water courses or drains.

14. TRANSPORT INFORMATION

Not classified as hazardous for transport.

15. REGULATORY INFORMATION

WHMIS: Not a controlled product

EPCRA 311/312 Categories: None

16. OTHER INFORMATION

This Material Safety Data Sheet (MSDS) may be used to comply with OSHA's Hazard Communication Standard, 29 CFR 1910.1200. Standard must be consulted for specific requirements.

To the best of our knowledge, the information contained herein is accurate. However, Roscoe Medical does not assume any liability whatsoever for the accuracy or completeness of the information contained herein. Final determination of suitability of any material is the sole responsibility of the user.

21973 Commerce Parkway
Strongsville, OH 44149

Toll free: 800-3-ROSCOE (376-7263)
main: 440-572-1962
fax: 440-572-42691

www.roscoemedical.com