



CIDEX[®] Tray Systems

INSTRUMENT HANDLING

CIDEX[®] Tray Systems have been designed by ASP to provide healthcare professionals with the ideal complementary accessories for manual cleaning/high-level disinfection. Manual processes have never been so easy.

EXPERIENCE THE BENEFITS

CIDEX[®] Tray Systems offer a large choice of premium durable trays.

- **Easy and practical handling of medical instruments**
- **Ideal complementary accessories for manual disinfection**
- **5 shapes and sizes to meet every user preference**



ASP[®] CIDEX[®] Tray Systems

FEATURES & BENEFITS

- **FLEXIBLE UTILIZATION** – can be used for either manual cleaning or high-level disinfection stages.
- **EASY PERIODIC STERILIZATION** – heat resistant up to 134°C.
- **LONG LIFE TIME** – superior durability of glass filled polypropylene.
- **REDUCTION OF VAPOR AND FUMES ESCAPE** with their unique rounded lip design.
- **EASY CLEANING** with their wide radius corner.
- **EASY SOLUTION MONITORING** – reusable labels help monitor the ASP High-level disinfection solutions reuse.

TECHNICAL INFORMATION

MEDICAL DEVICE CLASSIFICATION Class I according to MDD 93/42/EEC

ORDERING INFORMATION

REORDER NO.	DESCRIPTION AND DIMENSIONS	CONTENT
82010	344 mm x 177 mm x 132 mm 3.6 liters	1 tray, lid & liner
82016	503 mm x 186 mm x 132 mm 8.1 liters	1 tray, lid & liner
82027	740 mm x 220 mm x 132 mm 14.1 liters	1 tray, lid & liner
82032	449 mm diameter by 222 mm 17.5 liters	1 tray & lid
82076	594 mm x 381 mm x 141 mm 19.8 liters	1 tray, lid & liner



FOR MORE INFORMATION, CONTACT YOUR LOCAL ADVANCED STERILIZATION PRODUCTS SALES REPRESENTATIVE

www.aspjj.com/emea

<p>UK: Advanced Sterilization Products A Division of J&J Medical Ltd. Pinewood Campus Nine Mile Ride Wokingham Berkshire RG40, 3EW, England T: +44 1344 871 081 F: +44 1344 871 171</p>	<p>Ireland: Advanced Sterilization Products A Division of J&J Medical Ireland Airton Road, Tallaght, Dublin 24, Ireland T: +353 1 466 5200 F: +353 1 466 5340</p>	<p>Egypt: Johnson & Johnson Medical Egypt Florida Mall 5th Floor 1229 Square EL Sheikh Ali Gad EL Hak St. Heliopolis Cairo, Egypt T: +202 2268 5026 F: +202 2268 4674</p>	<p>Middle East: Advanced Sterilization Products A Division of J&J Medical Middle East Dubai Healthcare City, J&J building, 3rd floor PO Box 505080, Dubai, United Arab Emirates T: +971 4 429 7235 F: +971 4 429 7250</p>	<p>South Africa: Advanced Sterilization Products A Division of J&J Medical Pty. Ltd. SA PO Box 273, Midrand Halfway House 1685 South Africa T: +27 11 265 1120 F: +27 11 265 1189</p>
--	--	--	--	--

ADVANCED STERILIZATION PRODUCTS

Division of Medos International Sàrl
a **Johnson & Johnson** company