

## Clarity Urocheck 120 Urine Analyzer



- ✓ Up to 120 tests per hour
- ✓ Flagging for abnormal results
- ✓ 2000 Test Internal Memory
- ✓ Runs the Clarity 4OB, 5OB, URO7, and 10SG
- ✓ Automatic Paper Feeder
- ✓ Automatic Calibration and self check
- ✓ USB and RS232 Ports
- ✓ Barcode Reader (optional)
- ✓ Large LCD Display
- ✓ Competitively Priced

- ✓ Q.C. Lockout  
with Preprogrammed Values &  
Pass-Fail Summary
- ✓ Operator Lockout  
(Administrator & Operator  
Password Protection)
- ✓ Strip in same order as  
leading brand
- ✓ Fast & Easy To Use
- ✓ Patient ID Input  
via Manual or Bar Code Reader
- ✓ Compatible with Clarity's  
Middleware Solution  
for seamless data transfer to your  
EMR/LIS

**Product Codes: DTG-UAM, DTG-UAPROMO5\*, DTG-UAPROMO\*\***

*Contents: 1 Urocheck 120 Urine Analyzer. Operator Manual. Quick Start Guide.*

*\* Contains 5 Bottles of Urine Strips. \*\*Contains 10 Bottles of Urine Strips*

**Clear Vision - Clear Results**  
Clarity Brand Rapid Diagnostics

For more information call 877.722.6339 or 561.288.9028  
sales@claritydiagnostics.com      www.claritydiagnostics.com



## Clarity Urinalysis Strips

- ✓ Consistent
- ✓ Competitively Priced
- ✓ Strip in same order as leading brand



	DTG-10SG	DTG-URO7	DTG-5OB	DTG-4OB	DTG-2GP	DTG-KET	DTG-HEM
Leukocytes	X	X	X	X			
Nitrite	X	X	X				
Urobilinogen	X						
Protein	X	X	X	X	X		
pH	X	X					
Blood	X	X	X	X			X
Ketone	X	X				X	
Specific Gravity	X						
Bilirubin	X						
Glucose	X	X	X	X	X		

“The customer service is great. The quality of the **Clarity** strips are impeccable. The product works great with our existing analyzers at all of our urological centers.”

*Laboratory Director, Eyre Klindt - Metro Urology, St. Paul, MN*

Product Codes: DTG-10SG, DTG-URO7, DTG-5OB, DTG-4OB  
DTG-2GP, DTG-KET, DTG-HEM

Contents: 100 Urine Test Strips. 1 Package Insert.

**Clear Vision - Clear Results**  
Clarity Brand Rapid Diagnostics

For more information call 877.722.6339 or 561.288.9028  
sales@claritydiagnostics.com      www.claritydiagnostics.com

