

## Clarity Urinalysis Strips

- ✓ Consistent
- ✓ Competitively Priced
- ✓ Strip in same order as leading brand



	DTG-10SG	DTG-URO7	DTG-5OB	DTG-4OB	DTG-2GP	DTG-KET	CD-MAC25
Leukocytes	•	•	•	•			
Nitrites	•	•	•				
Urobilinogen	•						
Protein	•	•	•	•	•		
pH	•	•					
Blood	•	•	•	•			
Ketone	•	•				•	
Specific Gravity	•						
Bilirubin	•						
Glucose	•	•	•	•	•		
Microalbumin							•
Creatinine							•

Product Codes: DTG-10SG, DTG-URO7, DTG-5OB, DTG-4OB, DTG-2GP, DTG-KET, CD-MAC25\*

Contents: 100 Urine Test Strips. 1 Package Insert. \*CD-MAC25 contains 25 Urine Test Strips.

“The customer service is great. The quality of the **Clarity** strips are impeccable. The product works great with our existing analyzers at all of our urological centers.”

Laboratory Director, Eyre Klindt - Metro Urology, St. Paul, MN

**Clear Vision - Clear Results**  
Clarity Brand Rapid Diagnostics

For more information call 877.722.6339 or 561.288.9028  
sales@claritydiagnostics.com      www.claritydiagnostics.com

