



## Clarity Drug Testing Cups



- ✓ Integrated drug panel and testing cup in one.
- ✓ Simple one-step procedure.
- ✓ Screw lid for a no-leak cup.
- ✓ Easy-to-use.
- ✓ Generic product available.
- ✓ 510K cleared and CLIA Waived.

Drugs of Abuse	Abbreviation	Cutoff
Amphetamine	AMP	1000 ng/ml
Barbiturates	BAR	300 ng/ml
Buprenorphine/ Norbuprenorphine	BUP/NBUP	10 ng/ml
Benzodiazepine	BZO/BZD	300 ng/ml
Cocaine	COC	300 ng/ml
Methamphetamine	mAMP/MET	1000 ng/ml
MDMA or Ecstasy	MDMA	500 ng/ml
Morphine	MOP	300 ng/ml
Methadone	MTD	300 ng/ml
Opiate	OPI	2000 ng/ml
Oxycodone	OXY	100 ng/ml
Phencyclidine	PCP	25 ng/ml
Propoxyphene	PPX	300 ng/ml
Tricyclics	TCA	1000 ng/ml
Marijuana/Hashish	THC	50 ng/ml

Read the drug test results at 5 minutes, do not interpret the results after 10 minutes.



All positive results are presumptive and should be confirmed by an alternative method (e.g. GC/MS).


**Clear Vision - Clear Results**  
Clarity Brand Rapid Diagnostics

For more information call 877.722.6339 or 561.288.9028  
sales@claritydiagnostics.com www.claritydiagnostics.com





# Clarity Drug Testing Dip Cards




- ✓ Easy-to-use dip cards.
- ✓ Generic product available.
- ✓ 510K cleared and CLIA Waived.

Read the drug test results at 5 minutes, do not interpret the results after 10 minutes.



Drugs of Abuse	Abbreviation	Cutoff
Amphetamine	AMP	1000 ng/ml
Barbiturates	BAR	300 ng/ml
Buprenorphine/ Norbuprenorphine	BUP/NBUP	10 ng/ml
Benzodiazepine	BZO/BZD	300 ng/ml
Cocaine	COC	300 ng/ml
Methamphetamine	mAMP/MET	1000 ng/ml
MDMA or Ecstasy	MDMA	500 ng/ml
Morphine	MOP	300 ng/ml
Methadone	MTD	300 ng/ml
Opiate	OPI	2000 ng/ml
Oxycodone	OXY	100 ng/ml
Phencyclidine	PCP	25 ng/ml
Propoxyphene	PPX	300 ng/ml
Tricyclics	TCA	1000 ng/ml
Marijuana/Hashish	THC	50 ng/ml



All positive results are presumptive and should be confirmed by an alternative method (e.g. GC/MS).

**Clear Vision - Clear Results**  
Clarity Brand Rapid Diagnostics

For more information call 877.722.6339 or 561.288.9028  
sales@claritydiagnostics.com www.claritydiagnostics.com

